Effect of Audio-Visual Brain Entrainment on Mood and Quality of Sleep: a pilot trial with university students

Daiana Cristina Salm^{1,2}, Luiz Augusto Belmonte^{1,2}, Bruna Hoffmann de Oliveira^{1,2}, Luana Meneghini Belmonte^{1,2}, <u>Francisco J Cidral-Filho^{1,2}, Patrick Porter³ and Daniel Fernandes Martins^{1,2}</u>



¹ Experimental Neuroscience Laboratory (LaNEx), University of Southern Santa Catarina, Palhoça, Santa Catarina, Brazil. ² Postgraduate Program in Health Sciences, University of Southern Santa Catarina, Santa Catarina, Brazil. ³ Quantum University, Honolulu, HI - USA.

Contact: fcidral@braintap.com

OBJECTIVE

RESULTS

The Study Objective was to investigate the effect of the Audio-Visual Brain Entrainment

ABE effectively reduced ISI (data not statistically significant); PQSI (p<0.05); DASS-21



(ABE) on Mood and Quality of Sleep of university students.


METHODOLOGY

The study was conducted at the Laboratory of Experimental Neuroscience - University of Southern Santa Catarina (UNISUL), Brazil, and the protocol was approved by the Institutions Ethics committee. Informed consent forms were obtained during patient screening phase of the study at the site of the tests.


Sample size consisted of 7 university students (Four males and three females. Ages between 20 and 58), who were not making use of analgesics, anti-inflammatories or sleep aids 7 (seven) days prior to, as well as during the study, and who had no hearing disabilities.

ABE was delivered with a BrainTap headset (New Bern - NC – USA; Figure 1, A & B) in 20-minute sessions 3 times a week for 6 weeks. Session consists of Binaural beats at (data not statistically significant); and EPS-10 (data not statistically significant). The

participants reported feeling very relaxed during the sessions.



18 to 0.5 HZ, Isochronic Tones at 18 to .0.5 HZ and visual Entrainment through lightemitting diode lights at 470 nanometers (nm) flickering at 18 to 0.5 HZ. The following questionnaires were applied at baseline and after 6 weeks: Epworth Sleepiness Scale -Daytime sleepiness (ESS), Insomnia Severity Index (ISI), Pittsburgh Quality of Sleep Index (PQSI), Depression Anxiety and Stress Scale (DASS-21), and Perceived Stress Scale (EPS-10).






Figure 2 - Effect of Audio-Visual Brain Entrainment on Mood and Quality of Sleep. A) Epworth Sleepiness Scale; B) Insomnia Severity Index; C) Pittsburgh Quality of Sleep Index; D) Depression Anxiety and Stress Scale; E) Perceived Stress Scale. Data were expressed as mean \pm standard deviation (SD) n = 7. Student's T-test was used. *p<0.05 when compared with baseline evaluation. NS: not statistically significant.

Figure 1: A) Audio-visual Brain Entrainment Session. B) BrainTap headset.



ACKNOWLEDGMENT

Despite the reduced sample size (n=7), results indicate that Audio-Visual Brain

Entrainment (ABE) significantly increased Quality of Sleep of university students (PQSI

p<0.05). A lager sample size study is necessary to confirm and extend the effects of



2020 International Congress on **Integrative Medicine and Health** April 28 - May 1, 2020 | Cleveland, Ohio, USA Academic Consortium for Integrative Medicine and Health

Transforming Health Care Through Collaboration

ABE on Mood.